



ALFAPARF G R O U P

# YOUR SKIN IS REBORN WITH NEW LIGHT.

#### **PHOTOREJUVENATION**

Profound effect on fibroblasts with the production of new collagen and elastin at the dermis level.

#### SENSITIVE SKIN

Creates a selective thermal stimulus that performs discernible re-epithelialization of the blood vessel walls upon reaching them.

#### SKIN SPOTS

Significant reduction of the skin hyperchromias (hypermelanosis) caused by considerable increase of production of melanin by melanocytes.

### **ACNE AND GREASY SKIN**

Alleviates the inflammation process of the sebaceous gland, reducing the excessive production of sebum.





## TECHNICAL CHARACTERISTICS

INTENSE PULSED LIGHT EMITTER

DISPLAY: 81'CD touch screen DOUBLE EMISSION SYSTEM:

- single spot frequency at 1 Hz
- multi-frequency system: frequency from 1 Hz to 5 Hz TSF

(Top Speed Frequency)

HIGH DURABILITY XENON LAMP - 250,000 spots

EPILATION FILTERS: 3 filters (610 nm, 640 nm, 690 nm), plus

Double-Edged Light Emission Filter at 690-950 nm

MAXIMUM ENERGY DENSITY: 26 J/cm

PULSE DURATION: 2 ms - 20 ms

**SPOT AREA**: 50 mm X 15 mm = 7,5 cm

COOLING: with water using the double cooling system

(Peltier cells and copper radiators)