



DYNAMICWAVE
WITH **DOUBLE WAVE SYSTEM**

APG

TECH

Beauty in good hands.

ALFAPARF
GROUP

DYNAMIC WAVE WITH DOUBLE WAVE SYSTEM THE KO METHOD FOR CELLULITE*

DYNAMIC WAVE WITH DOUBLE WAVE SYSTEM

is a non-invasive technology that produces pressor acoustic waves that generate a compressing and decompressing action in the tissues, which is beneficial for fighting all stages of cellulite* and skin laxity.

DYNAMIC WAVE WITH DOUBLE WAVE SYSTEM

synergistically uses the energy produced by two specific handpieces :

- "BLUE WAVE" handpiece
- "RED WAVE" handpiece.

PLUS

- Increases the elasticity of the tissues.
- Firms the tissues.
- Redefines the silhouette.
- Enhances capillarity and vascularization.
- Oxygenates the tissues.
- Reduces the four stages of cellulite*.
- Improves the skin texture.
- Reduces the "orange-peel" appearance.

*Refers to the appearance of cellulite.



TECHNICAL CHARACTERISTICS

POWER SUPPLY: 230 VAC 50 Hz

MAXIMUM POWER CONSUMPTION: 700W

COMPRESSED AIR SUPPLY: 6 - 8 bar (600 - 800 kPa)
50 l/min at 8 bar (800 kPa)

OPERATING PRESSURE: 1,3 - 4 bar (130 - 400 kPa)

OPERATING MODE: multiple pulses with rate of repetition
8-12-16-20 Hz

DIAMETER OF THE TAPPING PLATES: 15mm, 36mm and 45mm

DIMENSIONS: Width 52cm, Depth 56cm, Height 115cm